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This instruction manual is only a supplement item for the HERMLE rotors operating manual, not a replacement! Please be sure to read the operating manual, before operating the rotors.


## Safety Notes

### ⚠ ATTENTION: Fixing the Rotor, Rotor Lids and Caps


- Before each run, please be sure that the rotor nut is tightened, rotor lids, caps and inserted carriers are correctly sealed.
- To secure the plastic lid with clip closure, it should first be inserted into the outer diameter of the rotor with one hand. Only then the plastic lid should be closed, by light pressure on the clip with the second hand. This guarantees a long life of the clip closure
- If unusual noise occurs upon the start of centrifugation; the rotor, the rotor lid or the cap may not be properly secured. Please stop centrifugation immediately, by pressing the STOP key.

### ⚠ ATTENTION: Loading Rotors

- Rotors must be loaded symmetrically and with equal weight.


Picture 1

Picture 2

- The adapter must only be loaded with the appropriate vessels.
- Always use the same type of vessels and plates, (material/density, volume and weight). As indicated by the manufacturer, always be aware of the max. g-force of the vessels and plates.
- For your safety, every bucket with the same rating number has to be placed on the opposite side of the rotor; this provides optimal performance and higher imbalance tolerance. The rating number is located on the top of each bucket.

Picture 3: Correct insertion of the buckets

Picture 4: Example of rating number on the bucket

**ATTENTION: Aggressive Chemicals**

- Rotors are high-quality components that withstand extremely heavy loads. To ensure the stability of the rotor, be aware of aggressive chemicals that can impair the stability of the rotor.
- Avoid the use of aggressive chemicals, including: strong weak alkali, strong acids, solutions with mercury, copper and other heavy metal ions, halogenated hydrocarbons, concentrated saline solutions and phenol.
- If the rotor, the rotor lids and /or the caps are contaminated with aggressive chemicals, clean it immediately, using a neutral cleaning agent and then rinse thoroughly with water. This particularly applies to the rotor bores.
- Rotor lids and caps, made of PC and PP, should be regularly checked for damage.
- Replace rotor lids or caps, which show cracks or milky stains, immediately.

**Lifetime**

Rotors and rotor lids, made of aluminum and stainless steel, have a maximum operating life of **7 years**, from first use. Transparent rotor lids and caps, made of PC and PP, as well as rotors, tube racks and adapters, made of PP, have a maximum operating life of up to **3 years**, from first use.

Conditions of the Operating Time:

-Proper use for a damage-free state, Recommended Care.

- Before each run, please check if the accessories are damaged. Replace any and all damaged accessories.
- Rotors, rotor lids, buckets, adapters or caps, which are showing any signs of corrosion or mechanical damage, are not longer functional.
- Do not use any accessories which are past their operating life.
- Be aware when inserting buckets and rotors, ensure that they do not become scratched.
- Protect the equipment from damage.

**Cleaning and Disinfection**

When cleaning rotors and accessories, (sleeves, adapters, buckets, rotor lids), use a mild cleaning agent. To prevent corrosion or material wear, we recommend cleaning at least once a week. Alongside with this, dirt/dust must be removed immediately.

1. Remove the rotor from the centrifuge.
2. Use a bottle brush to clean and disinfect the rotor bores.
3. Rotors and caps with sealing: Remove the seal and clean the groove.
4. Rinse the rotors, rotor lids and adapters with fresh water. Particularly, drilling of angle rotors.
5. For drying of the rotors and accessories, set them on a towel. Place the angle rotors with bores down to dry.
6. Dry the rotor cone with a soft, dry and lint-free cloth.
7. Do not grease the cone. Look for damage.
8. To prevent corrosion, due to moisture between the rotor and the motor shaft, the rotor should be removed at least once a week and the motor shaft should be cleaned with a clean, dry and lint-free cloth.

## Autoclaving

After Cleaning/Disinfecting, the rotor can be autoclaved, as described below.

The recommended time for autoclaving : 15-20 min at 121°C/249,8°F (2.15 bar).

- ⚠ ATTENTION: The sterilization time of 20 min. must not be exceeded! Continuous sterilization will cause a reduction of the mechanical resistance of the plastic material. We recommend to replace the seals of aerosol-tight lids and caps, after 5 autoclaving cycles.**

## Grease

1. If the swing out rotor does not swing freely after each time cleaning and autoclaving, apply a thin layer of rotor grease, (order number 38-5656), to the rotor bolts, (See **Picture 4**).
2. Ensure that the rotor cross pivots and the bucket grooves are free from contamination.

Picture 5


Picture 6

3. After each run of any fixed angle rotor, apply a thin layer of rotor grease, (Order Number 38-5656), to the lid thread, gasket and sealing surface, to prevent any damage, (see **Pictures 7 and 8**). This prevents the sealing from warping, during opening/closing. Additionally, ensure to grease the sealing of bucket caps resp. sealing surface regularly and before first use. (see **Pictures 6**). Otherwise, there is a risk that the seal dries out and thus the lid is difficult to open. To ensure tightness, be sure that the lid is tightly sealed.

Picture 7


Picture 8

## Replace the Seal

1. If necessary, use a blunt lever to lift the rubber seal from the nut. For this purpose, use e.g. the round side of a paper clip, (See **Picture 9**). Make sure that the seal will not become damaged.
2. Check the seal to make sure it is intact. Do not use any damaged, porous, discolored or dirty seals.
3. Insert the seal in one location along the groove, use your finger to spread it along the outer edge of the lid, until the seal is fully spread into the groove. Make sure that it's not waved.
4. Place the lid on the rotor and close the lid.
5. Remove the lid and check that the seal is positioned correctly.

**!** If the seal is too long or too short, remove the seal from the groove and re-insert it.  
**ATTENTION: Faulty sealing occurs when the gasket is handled incorrectly. Uniformly insert the seal. Do not pull the seal lengthwise!**

Picture 9

## Additional information for aerosol-tight rotors, lids and buckets

The aerosol tightness of rotors, rotor lids, buckets and caps has been tested and certified by the "TÜV Nord CERT GmbH, Certification Body Consumer Products, Essen (Germany)" in accordance with Annex AA of IEC 61010-2-020. The certificates can be downloaded on our webpage, [www.hermle-labortechnik.de](http://www.hermle-labortechnik.de). Aerosol-tight rotors and buckets are marked with the label "aerosol-tight".

**!** **ATTENTION: Autoclaving, mechanical stress and contamination by chemicals or other aggressive solvents can impair the aerosol-tightness of the rotors and buckets.**

- Check the integrity of the sealing and sealing surface of the aerosol-tight rotor lids or caps, before each use.
- Use only aerosol-tight rotor lids or caps if the seal is undamaged and clean.
- **Aerosol-tight plastic rotors as well as plastic lids and seals of aerosol-tight aluminum rotors and buckets must be replaced after 10 autoclaving cycles!**
- **Never** store closed aerosol-tight rotors or buckets.

HERMLE Labortechnik GmbH  
Siemensstraße 25  
78564 Wehingen  
Tel: +49-74 26-96 22-11  
Fax: +49-74 26-96 22-49  
Email: [info@hermlelt.de](mailto:info@hermlelt.de)  
Internet: <http://www.hermle-labortechnik.de>

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©HERMLE Labortechnik GmbH 2023