

Zubehör

Wird
ersetzt
in
Zentrifuge:

Z 216 M

Z 216 MK

Z 32 HK

Z 326

Z 326 K

Z 366

Z 366 K

Z 446

Z 446 K



NEU

Rotor 44 x 1.5/2.0 ml

Gerät	Z 216 M	Z 216 MK	Z 32 HK	Z 326	Z 326 K	Z 366	Z 366 K	Z 446	Z 446 K
ALTE Best. Nr.	220.88 V08	220.88 V08	220.88 V09	220.88 V09	220.88 V09	220.88 V09	220.88 V09	220.88 V09	220.88 V09
NEUE Best. Nr.	221.68 V01	221.68 V01	221.68 V02	221.68 V02	221.68 V02	221.68 V02	221.68 V02	221.68 V02	221.68 V02
max. Drehzahl	15000 min-1	15000 min-1	15000 min-1	15000 min-1	15000 min-1	15000 min-1	16000 min-1	16000 min-1	16000 min-1
max. RZB	21379 xg	21379 xg	21379 xg	21379 xg	21379 xg	21379 xg	24325 xg	24325 xg	24325 xg
Level Anlauf: 9/0 Auslauf: 9/0	35/259 18/188	32/265 18/258	25/221 22/476	28/257 21/174	28/256 21/434	31/256 24/327	33/275 26/536	25/236 25/324	25/235 25/500

704.004 Adapter Satz für 0.2/0.4 ml Röhrchen, Ø 6 mm

704.005 Adapter Satz für 0.5 ml Röhrchen, Ø 8 mm

HERMLE
LABORTECHNIK