

# Accessories

Now replaced in centrifuge:

**Z 216 M**

**Z 216 MK**

**Z 32 HK**

**Z 326**

**Z 326 K**

**Z 366**

**Z 366 K**

**Z 446**

**Z 446 K**



# NEW

## Rotor 44 x 1.5/2.0 ml

| Unit                 | Z 216 M           | Z 216 MK          | Z 32 HK           | Z 326             | Z 326 K           | Z 366             | Z 366 K           | Z 446             | Z 446 K           |
|----------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|
| OLD Order No.        | 220.88 V08        | 220.88 V08        | 220.88 V09        | 220.88 V09        | 220.88 V09        | 220.88 V09        | 220.88 V09        | 220.88 V09        | 220.88 V09        |
| <b>NEW Order No.</b> | <b>221.68 V01</b> | <b>221.68 V01</b> | <b>221.68 V02</b> | <b>221.68 V02</b> | <b>221.68 V02</b> | <b>221.68 V02</b> | <b>221.68 V02</b> | <b>221.68 V02</b> | <b>221.68 V02</b> |
| max. Speed           | 15000 rpm         | 15000 rpm         | 15000 rpm         | 15000 rpm         | 15000 rpm         | 15000 rpm         | 16000 rpm         | 16000 rpm         | 16000 rpm         |
| max. RCF             | 21379 xg          | 21379 xg          | 21379 xg          | 21379 xg          | 21379 xg          | 21379 xg          | 24325 xg          | 24325 xg          | 24325 xg          |
| Level                |                   |                   |                   |                   |                   |                   |                   |                   |                   |
| Acc: 9/0             | 35/259            | 32/265            | 25/221            | 28/257            | 28/256            | 31/256            | 33/275            | 25/236            | 25/235            |
| Dec: 9/0             | 18/188            | 18/258            | 22/476            | 21/174            | 21/434            | 24/327            | 26/536            | 25/324            | 25/500            |

**704.004** Adapter set for 0.2/0.4 ml tubes, Ø 6 mm

**704.005** Adapter set for 0.5 ml tubes, Ø 8 mm